FRENCH EDITION

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
Courbe de survie de Kaplan-Meier

References

1. **Reigstad O, Holm-Glad T, Bolstad B, Grimsgaard C, Thorkildsen R, Røkkum M.** Five- to 10-Year Prospective Follow-Up of Wrist Arthroplasty in 56 Nonrheumatoid Patients. *J Hand Surg Am.* 2017
2. **Glad TH, Røkkum M, Thorkildsen R, Reigstad O.** Rearticulation from arthrodesis with Motec prosthesis (Reartikulering av avstivede håndledd med Motec protese). *NOF Høstmøteboken Abstract 210*, 2016.
3. **Youm Y, McMurthy RY, Flatt AE, Gillespie TE.** Kinematics of the wrist. I. An experimental study of radial-ulnar deviation and flexion-extension. *J Bone Joint Surg Am.* 1978 Jun;60(4):423-31.
4. **Utzschneider S, Becker F, Grupp TM, Sievers B, Paulus A, Gottschalk O, Jansson V.** Inflammatory response against different carbon fiber-reinforced PEEK wear particles compared with UHMWPE in vivo. *Acta Biomater.* 2010 Nov;6(11):4296-304.
5. **Reigstad A, Reigstad O, Grimsgaard C, Røkkum M.** New concept for total wrist replacement. *J Plast Surg Hand Surg.* 2011 Jun;45(3):148-56.
6. **Glad TH, Røkkum M, Thorkildsen R, Reigstad O.** Chrome and Cobalt blood levels in a metal-on-metal total wrist arthroplasty. *NOF Høstmøteboken Abstract 212*, 2016.
7. **Lundborg G, Besjakov J, Brånemark PI.** Osseointegrated wrist-joint prostheses: a 15-year follow-up with focus on bony fixation. *Scand J Plast Reconstr Surg Hand Surg.* 2007;41(3):130-7.
8. **Røkkum M.** On Late Complications With Ha Coated Hip Arthroplasties, Department of Biomaterials/Handicap Research, Institute for Surgical Sciences, Faculty of Medicine, University of Göteborg, Göteborg, Sweden and Orthopaedic University Clinic, National Hospital, Oslo, Norway, Göteborg 2001.
9. **Reigstad O, Franke-Stenport V, Johansson CB, Wennerberg A, Røkkum M, Reigstad A.** Improved bone ingrowth and fixation with a thin calcium phosphate coating intended for complete resorption. *J Biomed Mater Res B Appl Biomater.* 2007 Oct;83(1):9-15.
10. **Reigstad O, Johansson C, Stenport V, Wennerberg A, Reigstad A, Røkkum M.** Different patterns of bone fixation with hydroxyapatite and resorbable CaP coatings in the rabbit tibia at 6, 12, and 52 weeks. *J Biomed Mater Res B Appl Biomater.* 2011 Oct;99(1):14-20.
11. **Reigstad O, Lütken T, Grimsgaard C, Bolstad B, Thorkildsen R, Røkkum M.** Promising one- to six-year results with the Motec wrist arthroplasty in patients with post-traumatic osteoarthritis. *J Bone Joint Surg Br.* 2012 Nov;94(11):1540-5.
12. **Reigstad O, Røkkum M.** Conversion of Total Wrist Arthroplasty to Arthrodesis with a Custom-Made Peg *J Wrist Surg* 2014;3:211–215
13. **Scholes SC, Unsworth A.** Pitch-based carbon-fibre-reinforced poly (ether-ether-ketone) OPTIMA assessed as a bearing material in a mobile bearing unicondylar knee joint. *Proc Inst Mech Eng H.* 2009 Jan;223(1):13-25.
14. **Scholes SC, Unsworth A.** Wear studies on the likely performance of CFR-PEEK/CoCrMo for use as artificial joint bearing materials. *J Mater Sci Mater Med.* 2009 Jan;20(1):163-70.
15. **Kabir K, Schwiesau J, Burger C, Pflugmacher R, Grupp T, Wirtz DC.** Comparison of Biological Response to UHMWPE and CFR-PEEK Particles in Epidural Space. *Universitätsklinikum Bonn, Department for Orthopaedics and Trauma Surgery, Bonn, Germany, Aesculap AG Research & Development, Tuttlingen, Germany.* 2011.
16. **Grupp TM, Utzschneider S, Schröder C, Schwiesau J, Fritz B, Maas A, Blömer W, Jansson V.** Biotribology of alternative bearing materials for unicompartamental knee arthroplasty. *Acta Biomater.* 2010 Sep;6(9):3601-10. Erratum in: *Acta Biomater.* 2012 Apr;8(4):1659.
17. **Wennerberg A.** On surface roughness and implant incorporation. *Department of Biomaterial/Handicap Research, Göteborg, Sweden.* 1996.

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


Product information






- Needed for CFR-PEEK articulation
- Needed for CoCrMo articulation

Implants













Radius Cup CoCrMo Ø15 mm	40-1015S	●	
Radius Cup CoCrMo Ø18 mm (optional)	40-1018S	●	
Radius Cup CFR-PEEK Ø15 mm	40-1915S	●	
Metacarpal Head Ø15 mm Short Neck (optional)	40-1115S	● ●	
Metacarpal Head Ø18 mm Short Neck (optional)	40-1118S	●	
Metacarpal Head Ø15 mm Medium Neck	40-1715S	● ●	
Metacarpal Head Ø18 mm Medium Neck (optional)	40-1718S	●	
Metacarpal Head Ø15 mm Long Neck	40-1215S	● ●	
Metacarpal Head Ø18 mm Long Neck (optional)	40-1218S	●	
Metacarpal Head Ø15 mm Extra Long Neck	40-1315S	● ●	
Radius Threaded Implant length 32 mm	40-1332S	● ●	
Radius Threaded Implant length 38 mm	40-1338S	● ●	
Radius Threaded Implant length 44 mm	40-1344S	● ●	
Radius Threaded Implant length 50 mm	40-1350S	● ●	
Radius Threaded Implant length 56-80 mm (optional)			
Metacarpal Threaded Implant length 45 mm Large	40-1445S	● ●	
Metacarpal Threaded Implant length 50 mm Large	40-1450S	● ●	
Metacarpal Threaded Implant length 55 mm Large	40-1455S	● ●	
Metacarpal Threaded Implant length 60 mm Large	40-1460S	● ●	

Metacarpal Threaded Implant length 65 mm Large (optional)	40-1465S	 	
Metacarpal Threaded Implant length 70 mm Large (optional)	40-1470S	 	
Metacarpal Threaded Implant length 45 mm Small	40-1475S	 	
Metacarpal Threaded Implant length 50 mm Small	40-1480S	 	
Metacarpal Threaded Implant length 55 mm Small	40-1485S	 	
Metacarpal Threaded Implant length 60 mm Small	40-1490S	 	
Metacarpal Threaded Implant length 65 mm Small (optional)	40-1495S	 	
Metacarpal Threaded Implant length 70 mm Small (optional)	40-1400S	 	

Trials

Trial – Radius Cup Ø15 mm	40-1522		
Trial – Radius Cup Ø18 mm (optional)	40-1521		
Trial – Radius Cup Ø15 mm For CFR-PEEK Cup	40-1541		
Trial – Metacarpal Head Ø15 mm Short Neck (optional)	40-1529	 	
Trial – Metacarpal Head Ø18 mm Short Neck (optional)	40-1527		
Trial – Metacarpal Head Ø15 mm Medium Neck	40-1524	 	
Trial – Metacarpal Head Ø18 mm Medium Neck (optional)	40-1523		
Trial – Metacarpal Head Ø15 mm Long Neck	40-1528	 	
Trial – Metacarpal Head Ø18 mm Long Neck (optional)	40-1526		
Trial – Metacarpal Head Ø15 mm Extra Long Neck	40-1602	 	

Instruments

Hohmann Retractor	40-1503	● ●	
Bits 3,5 mm HEX with Quick-Lock	40-1513	● ●	
Impactor	40-1516	● ●	
Guide Wire T-handle	40-1518	● ●	
Cup Remover	40-1519	● ●	
Cannulated Drill for Radius 32-50 mm	40-1546	● ●	
Cannulated Drill for Metacarpal III 45-70 mm Large	40-1551	● ●	
Cannulated Drill for Metacarpal III 45-70 mm Small	40-1552	● ●	
Guide Wire with sharp tip Ø2 mm	40-1561	● ●	
Guide Wire with round tip Ø2 mm	40-1563	● ●	
Radius Spherical Drill Ø18 mm (optional)	40-1566	●	
Radius Spherical Drill Ø15 mm	40-1567	● ●	
Handle Tri-Lobe with Quick-Lock	45-2585	● ●	
Handle Tri-Lobe with Ratchet (optional)	40-2593	● ●	
Adapter, from AO male to Tri-Lobe female (optional)	40-5000	● ●	
Hammer	52-2211	● ●	
Awl	62-3070	● ●	
Tray and lid	40-1600	● ●	

IFU

For the latest version of this Instruction For Use. Please visit:
<http://download.swemac.com/Motec-Wrist-Joint-Prosthesis>

Swemac develops and promotes innovative solutions for fracture treatment and joint replacement. We create outstanding value for our clients and their patients by being a very competent and reliable partner.

Swemac

Motec Wrist Joint Prosthesis

Manufacturer

Swemac Innovation AB

CE 0413

Cobolgatan 1 • SE-583 35 Linköping • Sweden

+46 13 37 40 30 • info@swemac.com • www.swemac.com

P125-28-1-FR-20190227

This brochure is a translation of: P125-28-1-20190226