GERMAN EDITION

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
[1] Costi, J. und Bain, G. Biomechanical Assessment of the Motec Wrist Joint Prosthesis, Report ID: 101218-BIL, 2018.

[2] Fraysse, F. et al., 2014. A novel method to replicate the kinematics of the carpus using a six degree-of-freedom robot. Journal of Biomechanics.

[3] Scott W. Wolfe et al. The Dart-Throwing Motion of the Wrist: Is it Unique to Humans? J Hand Surg Am, 2006.

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
Reigstad O., Røkkum M.:

Conversion of Total Wrist Arthroplasty to Arthrodesis with a Custom-Made Peg
J Wrist Surg 2014; 3; 211–215

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



Five- to 10-Year Prospective Follow-Up of Wrist Arthroplasty in 56 Nonrheumatoid Patients

Reigstad O, Holm-Glad T, Bolstad B, Grimsgaard C, Thorkildsen R, Røkkum M.


J Hand Surg Am. 2017. Oct;42(10):788-796.

Abstract

Fifty-seven (40 male) patients with end-stage arthritis changes received an uncemented ball-and-socket total wrist arthroplasty (Motec Wrist). Function was evaluated before surgery and at yearly follow-ups. Visual analog scale at rest and activity, quick Disabilities of the Arm, Shoulder, and Hand (QuickDASH), active range of motion (AROM), and grip-strength were recorded. Standardized radiographs were taken to assess osteolysis, loosening, and subsidence.

Fifty-six patients were followed for a mean of 8 years (SD, 2 years). Eight wrists were reoperated with arthrodesis (4) or a new arthroplasty (4) owing to distal component loosening (3), infection (2), pain/fixed malposition (2), or proximal and distal component loosening (1). One radiocarpal dislocation was reduced closed and remained stable. Improved QuickDASH score and visual analog scale pain score both at rest and during activity were found at the last follow-up, as well as increased AROM (97 vs 126) and grip strength (21 kg vs 24 kg). The radiological follow-up demonstrated loosening in 2 wrists. Thirty-five patients were working at surgery (17 manual labor) and 27 (11 manual labor) at follow-up. The 10-year Kaplan-Meier survival of the implants was 86% for revision any cause, 2 additional arthroplasties are loose (but not revised), giving a survival rate of 82% if these are revised prior to 10 years of observation.

An uncemented total wrist arthroplasty can provide long-lasting unrestricted hand function in young and active patients.

Kaplan-Meier survival curve








Literatur

1. **Reigstad O., Røkkum M.:** Wrist arthroplasty using prosthesis as an alternative to arthrodesis: design, outcomes and future. *J Hand Surg Eur Vol.* Sep 2018; 43(7); 689-699. doi: 10.1177/1753193418784707
2. **Reigstad O., Holm-Glad T., Bolstad B., Grimsgaard C., Thorkildsen R., Røkkum M.:** Five- to 10-Year Prospective Follow-Up of Wrist Arthroplasty in 56 Nonrheumatoid Patients. *J Hand Surg Am.* Okt 2017; 42(10); 788-796.
3. **Reigstad A., Reigstad O., Grimsgaard C., Røkkum M.:** New concept for total wrist replacement. *J Plast Surg Hand Surg.* Juni 2011; 45(3); 148-56.
4. **Reigstad O., Lütken T., Grimsgaard C., Bolstad B., Thorkildsen R., Røkkum M.:** Promising one- to six-year results with the Motec wrist arthroplasty in patients with post-traumatic osteoarthritis. *J Bone Joint Surg Br.* Nov 2012; 94(11); 1540-5.
5. **Holm-Glad T., Røkkum M., Thorkildsen R., Reigstad O.:** Rearticulation from arthrodesis with Motec prosthesis (Reartikulering av avstivede håndledd med Motec protese). *NOF Høstmøteboken Abstract 210*, 2016.
6. **Youm Y., McMurthy R.Y., Flatt A.E., Gillespie T.E.:** Kinematics of the wrist. I. An experimental study of radial-ulnar deviation and flexion-extension. *J Bone Joint Surg Am.* Jun 1978; 60(4); 423-31.
7. **Utzschneider S., Becker F., Grupp T.M., Sievers B., Paulus A., Gottschalk O., Jansson V.:** Inflammatory response against different carbon fiber-reinforced PEEK wear particles compared with UHMWPE in vivo. *Acta Biomater.* Nov 2010; 6(11); 4296-304.
8. **Scholes S.C., Unsworth A.:** Pitch-based carbon-fibre-reinforced poly (ether-ether-ketone) OPTIMA assessed as a bearing material in a mobile bearing unicondylar knee joint. *Proc Inst Mech Eng H.* Jan 2009; 223(1); 13-25.
9. **Scholes S.C., Unsworth A.:** Wear studies on the likely performance of CFR-PEEK/CoCrMo for use as artificial joint bearing materials. *J Mater Sci Mater Med.* Jan 2009; 20(1); 163-70.
10. **Kabir K., Schwiesau J., Burger C., Pflugmacher R., Grupp T., Wirtz D.C.:** Comparison of Biological Response to UHMWPE and CFR-PEEK Particles in Epidural Space. *Universitätsklinikum Bonn, Orthopädie und Unfallchirurgie, Bonn, Deutschland, Aesculap AG Research & Development, Tuttlingen, Deutschland.* 2011.
11. **Grupp T.M., Utzschneider S., Schröder C., Schwiesau J., Fritz B., Maas A., Blömer W., Jansson V.:** Biotribology of alternative bearing materials for unicompartamental knee arthroplasty. *Acta Biomater.* Sep 2010; 6(9); 3601-10. Erratum in: *Acta Biomater.* Apr 2012; 8(4); 1659.
12. **Holm-Glad T., Røkkum M., Thorkildsen R., Grimsgaard C., Reigstad O.:** Low chrome and cobalt blood levels with a metal-on-metal total wrist arthroplasty. *FESSH 2014.* Hand- and microsurgery section, Orthopaedic department, OUS-Rikshospitalet, Oslo, Norwegen.
13. **Reigstad O.:** Wrist arthroplasty- Bone fixation, clinical development and mid to long term, *Acta Orthopaedica (Suppl 354)* 2014; 85.
14. **Lundborg G., Besjakov J., Brånemark P.I.:** Osseointegrated wrist-joint prostheses: a 15-year follow-up with focus on bony fixation. *Scand J Plast Reconstr Surg Hand Surg.* 2007; 41(3); 130-7.
15. **Wennerberg A.:** On surface roughness and implant incorporation. *Department of Biomaterial/Handicap Research, Göteborg, Schweden.* 1996.
16. **Reigstad O., Franke-Stenport V., Johansson C.B., Wennerberg A., Røkkum M., Reigstad A.:** Improved bone ingrowth and fixation with a thin calcium phosphate coating intended for complete resorption. *J Biomed Mater Res B Appl Biomater.* Okt 2007; 83(1); 9-15.
17. **Reigstad O., Johansson C., Stenport V., Wennerberg A., Reigstad A., Røkkum M.:** Different patterns of bone fixation with hydroxyapatite and resorbable CaP coatings in the rabbit tibia at 6, 12, and 52 weeks. *J Biomed Mater Res B Appl Biomater.* Okt 2011; 99(1); 14-20.
18. **Costi J., Bain G.:** Biomechanical Assessment of the Motec Wrist Joint Prosthesis, Report ID: 101218-BIL, 2018.
19. **Reigstad O., Røkkum M.:** Conversion of Total Wrist Arthroplasty to Arthrodesis with a Custom-Made Peg *J Wrist Surg* 2014; 3; 211–215.

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











Metacarpal III Threaded Implant length 65 mm Large	40-1465S	● ●	
Metacarpal III Threaded Implant length 70 mm Large	40-1470S	● ●	
Metacarpal III Threaded Implant length 45 mm Small	40-1475S	● ●	
Metacarpal III Threaded Implant length 50 mm Small	40-1480S	● ●	
Metacarpal III Threaded Implant length 55 mm Small	40-1485S	● ●	
Metacarpal III Threaded Implant length 60 mm Small	40-1490S	● ●	
Metacarpal III Threaded Implant length 65 mm Small	40-1495S	● ●	
Metacarpal III Threaded Implant length 70 mm Small	40-1400S	● ●	

Probeimplantate

Trial – Radius Cup Ø15 mm	40-1522	●	
Trial – Radius Cup Ø18 mm	40-1521	●	
Trial – Radius Cup Ø15 mm For CFR-PEEK Cup	40-1541	●	
Trial – Metacarpal Head Ø15 mm Short Neck	40-1529	● ●	
Trial – Metacarpal Head Ø18 mm Short Neck	40-1527	●	
Trial – Metacarpal Head Ø15 mm Medium Neck	40-1524	● ●	
Trial – Metacarpal Head Ø18 mm Medium Neck	40-1523	●	
Trial – Metacarpal Head Ø15 mm Long Neck	40-1528	● ●	
Trial – Metacarpal Head Ø18 mm Long Neck	40-1526	●	
Trial – Metacarpal Head Ø15 mm Extra Long Neck	40-1602	● ●	

Instrumente

Hohmann Capitate Retractor	40-1503	● ●	
Bits 3,5 mm HEX with Quick-Lock	40-1513	● ●	
Impactor	40-1516	● ●	
Guide Wire T-handle	40-1518	● ●	
Cup Remover	40-1519	● ●	
Cannulated Drill for Radius 32-50 mm	40-1546	● ●	
Cannulated Drill for Metacarpal III 45-70 mm Large	40-1551	● ●	
Cannulated Drill for Metacarpal III 45-70 mm Small	40-1552	● ●	
Guide Wire with sharp tip Ø2 mm	40-1561	● ●	
Guide Wire with round tip Ø2 mm	40-1563	● ●	
Radius Spherical Drill Ø18 mm	40-1566	●	
Radius Spherical Drill Ø15 mm	40-1567	● ●	
Handle Tri-Lobe with Quick-Lock	45-2585	● ●	
Handle Tri-Lobe with Ratchet (optional)	40-2593	● ●	
Adapter, from AO male to Tri-Lobe female (optional)	40-5000	● ●	
Hammer	52-2211	● ●	
Awl	62-3070	● ●	
Tray and lid	40-1600	● ●	

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Swemac develops and promotes innovative solutions for fracture treatment and joint replacement. We create outstanding value for our clients and their patients by being a very competent and reliable partner.

Swemac

Motec Wrist Joint Prosthesis

Hersteller

Swemac Innovation AB

CE 0413

Cobolgatan 1 • SE-583 35 Linköping • Schweden

+46 13 37 40 30 • info@swemac.com • www.swemac.com

P125-28-1-DE-20200310

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