FRENCH EDITION

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
Reigstad O, Røkkum M.

Conversion of Total Wrist Arthroplasty to Arthrodesis with a Custom-Made Peg
J Wrist Surg 2014;3:211-215

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


Five- to 10-Year Prospective Follow-Up of Wrist Arthroplasty in 56 Nonrheumatoid Patients

Reigstad O, Holm-Glad T, Bolstad B, Grimsgaard C, Thorkildsen R, Røkkum M.


J Hand Surg Am. 2017. Oct;42(10):788-796.

Abstract

Fifty-seven (40 male) patients with end-stage arthritis changes received an uncemented ball-and-socket total wrist arthroplasty (Motec Wrist). Function was evaluated before surgery and at yearly follow-ups. Visual analog scale at rest and activity, quick Disabilities of the Arm, Shoulder, and Hand (QuickDASH), active range of motion (AROM), and grip-strength were recorded. Standardized radiographs were taken to assess osteolysis, loosening, and subsidence.

Fifty-six patients were followed for a mean of 8 years (SD, 2 years). Eight wrists were reoperated with arthrodesis (4) or a new arthroplasty (4) owing to distal component loosening (3), infection (2), pain/fixed malposition (2), or proximal and distal component loosening (1). One radiocarpal dislocation was reduced closed and remained stable. Improved QuickDASH score and visual analog scale pain score both at rest and during activity were found at the last follow-up, as well as increased AROM (97 vs 126) and grip strength (21 kg vs 24 kg). The radiological follow-up demonstrated loosening in 2 wrists. Thirty-five patients were working at surgery (17 manual labor) and 27 (11 manual labor) at follow-up. The 10-year Kaplan-Meier survival of the implants was 86% for revision any cause, 2 additional arthroplasties are loose (but not revised), giving a survival rate of 82% if these are revised prior to 10 years of observation.

An uncemented total wrist arthroplasty can provide long-lasting unrestricted hand function in young and active patients.

Kaplan-Meier survival curve

Références




1. **Reigstad O, Røkkum M.** Wrist arthroplasty using prosthesis as an alternative to arthrodesis: design, outcomes and future. *J Hand Surg Eur Vol.* 2018 Sep;43(7):689-699. doi: 10.1177/1753193418784707
2. **Reigstad O, Holm-Glad T, Bolstad B, Grimsgaard C, Thorkildsen R, Røkkum M.** Five- to 10-Year Prospective Follow-Up of Wrist Arthroplasty in 56 Nonrheumatoid Patients. *J Hand Surg Am.* 2017. Oct;42(10):788-796.
3. **Reigstad A, Reigstad O, Grimsgaard C, Røkkum M.** New concept for total wrist replacement. *J Plast Surg Hand Surg.* 2011 Jun;45(3):148-56.
4. **Reigstad O, Lütken T, Grimsgaard C, Bolstad B, Thorkildsen R, Røkkum M.** Promising one- to six-year results with the Motec wrist arthroplasty in patients with post-traumatic osteoarthritis. *J Bone Joint Surg Br.* 2012 Nov;94(11):1540-5.
5. **Holm-Glad T, Røkkum M, Thorkildsen R, Reigstad O.** Rearticulation from arthrodesis with Motec prosthesis (Reartikulering av avstivede håndledd med Motec protese). *NOF Høstmøteboken Abstract 210*, 2016.
6. **Youm Y, McMurthy RY, Flatt AE, Gillespie TE.** Kinematics of the wrist. I. An experimental study of radial-ulnar deviation and flexion-extension. *J Bone Joint Surg Am.* 1978 Jun;60(4):423-31.
7. **Utzschneider S, Becker F, Grupp TM, Sievers B, Paulus A, Gottschalk O, Jansson V.** Inflammatory response against different carbon fiber-reinforced PEEK wear particles compared with UHMWPE in vivo. *Acta Biomater.* 2010 Nov;6(11):4296-304.
8. **Scholes SC, Unsworth A.** Pitch-based carbon-fibre-reinforced poly (ether-ether-ketone) OPTIMA assessed as a bearing material in a mobile bearing unicondylar knee joint. *Proc Inst Mech Eng H.* 2009 Jan;223(1):13-25.
9. **Scholes SC, Unsworth A.** Wear studies on the likely performance of CFR-PEEK/CoCrMo for use as artificial joint bearing materials. *J Mater Sci Mater Med.* 2009 Jan;20(1):163-70.
10. **Kabir K, Schwiesau J, Burger C, Pflugmacher R, Grupp T, Wirtz DC.** Comparison of Biological Response to UHMWPE and CFR-PEEK Particles in Epidural Space. *Universitätsklinikum Bonn, Department for Orthopaedics and Trauma Surgery, Bonn, Germany, Aesculap AG Research & Development, Tuttlingen, Germany.* 2011.
11. **Grupp TM, Utzschneider S, Schröder C, Schwiesau J, Fritz B, Maas A, Blömer W, Jansson V.** Biotribology of alternative bearing materials for unicompartmental knee arthroplasty. *Acta Biomater.* 2010 Sep;6(9):3601-10. Erratum in : *Acta Biomater.* 2012 Apr;8(4):1659.
12. **Holm-Glad T, Røkkum M, Thorkildsen R, Grimsgaard C and Reigstad O.** Low chrome and cobalt blood levels with a metal-on-metal total wrist arthroplasty. *FESSH 2014.* Hand- and microsurgery section, Orthopaedic department, OUS-Rikshospitalet, Oslo, Norway.
13. **Reigstad O.** Wrist arthroplasty- Bone fixation, clinical development and mid to long term, *Acta Orthopaedica (Suppl 354)* 2014; 85.
14. **Lundborg G, Besjakov J, Brånemark PI.** Osseointegrated wrist-joint prostheses: a 15-year follow-up with focus on bony fixation. *Scand J Plast Reconstr Surg Hand Surg.* 2007;41(3):130-7.
15. **Wennerberg A.** On surface roughness and implant incorporation. *Department of Biomaterial/Handicap Research, Göteborg, Sweden.* 1996.
16. **Reigstad O, Franke-Stenport V, Johansson CB, Wennerberg A, Røkkum M, Reigstad A.** Improved bone ingrowth and fixation with a thin calcium phosphate coating intended for complete resorption. *J Biomed Mater Res B Appl Biomater.* 2007 Oct;83(1):9-15.
17. **Reigstad O, Johansson C, Stenport V, Wennerberg A, Reigstad A, Røkkum M.** Different patterns of bone fixation with hydroxyapatite and resorbable CaP coatings in the rabbit tibia at 6, 12, and 52 weeks. *J Biomed Mater Res B Appl Biomater.* 2011 Oct;99(1):14-20.
18. **Costi, J. and Bain, G.** Biomechanical Assessment of the Motec Wrist Joint Prosthesis, Report ID: 101218-BIL, 2018.
19. **Reigstad O, Røkkum M.** Conversion of Total Wrist Arthroplasty to Arthrodesis with a Custom-Made Peg *J Wrist Surg* 2014;3:211–215

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
● Requis pour articulation CFR-PEEK







● Requis pour articulation CoCrMo

Implants









Radius Cup CoCrMo Ø15 mm	40-1015S	●	
Radius Cup CoCrMo Ø18 mm	40-1018S	●	
Radius Cup CFR-PEEK Ø15 mm	40-1915S	●	
Metacarpal Head Ø15 mm Short Neck	40-1115S	● ●	
Metacarpal Head Ø18 mm Short Neck	40-1118S	●	
Metacarpal Head Ø15 mm Medium Neck	40-1715S	● ●	
Metacarpal Head Ø18 mm Medium Neck	40-1718S	●	
Metacarpal Head Ø15 mm Long Neck	40-1215S	● ●	
Metacarpal Head Ø18 mm Long Neck	40-1218S	●	
Metacarpal Head Ø15 mm Extra Long Neck	40-1315S	● ●	
Radius Threaded Implant length 32 mm	40-1332S	● ●	
Radius Threaded Implant length 38 mm	40-1338S	● ●	
Radius Threaded Implant length 44 mm	40-1344S	● ●	
Radius Threaded Implant length 50 mm	40-1350S	● ●	
Metacarpal III Threaded Implant length 45 mm Large	40-1445S	● ●	
Metacarpal III Threaded Implant length 50 mm Large	40-1450S	● ●	
Metacarpal III Threaded Implant length 55 mm Large	40-1455S	● ●	
Metacarpal III Threaded Implant length 60 mm Large	40-1460S	● ●	

Metacarpal III Threaded Implant length 65 mm Large	40-1465S	● ●	
Metacarpal III Threaded Implant length 70 mm Large	40-1470S	● ●	
Metacarpal III Threaded Implant length 45 mm Small	40-1475S	● ●	
Metacarpal III Threaded Implant length 50 mm Small	40-1480S	● ●	
Metacarpal III Threaded Implant length 55 mm Small	40-1485S	● ●	
Metacarpal III Threaded Implant length 60 mm Small	40-1490S	● ●	
Metacarpal III Threaded Implant length 65 mm Small	40-1495S	● ●	
Metacarpal III Threaded Implant length 70 mm Small	40-1400S	● ●	

Implants d'essai

Trial – Radius Cup Ø15 mm	40-1522	●	
Trial – Radius Cup Ø18 mm	40-1521	●	
Trial – Radius Cup Ø15 mm For CFR-PEEK Cup	40-1541	●	
Trial – Metacarpal Head Ø15 mm Short Neck	40-1529	● ●	
Trial – Metacarpal Head Ø18 mm Short Neck	40-1527	●	
Trial – Metacarpal Head Ø15 mm Medium Neck	40-1524	● ●	
Trial – Metacarpal Head Ø18 mm Medium Neck	40-1523	●	
Trial – Metacarpal Head Ø15 mm Long Neck	40-1528	● ●	
Trial – Metacarpal Head Ø18 mm Long Neck	40-1526	●	
Trial – Metacarpal Head Ø15 mm Extra Long Neck	40-1602	● ●	

Instruments

Hohmann Capitate Retractor	40-1503	● ●	
Bits 3,5 mm HEX with Quick-Lock	40-1513	● ●	
Impactor	40-1516	● ●	
Guide Wire T-handle	40-1518	● ●	
Cup Remover	40-1519	● ●	
Cannulated Drill for Radius 32-50 mm	40-1546	● ●	
Cannulated Drill for Metacarpal III 45-70 mm Large	40-1551	● ●	
Cannulated Drill for Metacarpal III 45-70 mm Small	40-1552	● ●	
Guide Wire with sharp tip Ø2 mm	40-1561	● ●	
Guide Wire with round tip Ø2 mm	40-1563	● ●	
Radius Spherical Drill Ø18 mm	40-1566	● ●	
Radius Spherical Drill Ø15 mm	40-1567	● ●	
Handle Tri-Lobe with Quick-Lock	45-2585	● ●	
Handle Tri-Lobe with Ratchet (optional)	40-2593	● ●	
Adapter, from AO male to Tri-Lobe female (optional)	40-5000	● ●	
Hammer	52-2211	● ●	
Awl	62-3070	● ●	
Tray and lid	40-1600	● ●	

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Swemac develops and promotes innovative solutions for fracture treatment and joint replacement. We create outstanding value for our clients and their patients by being a very competent and reliable partner.

Swemac

Motec Wrist Joint Prosthesis

Fabricant

Swemac Innovation AB

CE 0413

Cobolgatan 1 • SE-583 35 Linköping • Suède

+46 13 37 40 30 • info@swemac.com • www.swemac.com

P125-28-1-FR-20200303

Traduction de : P125-28-1-20191204