

TOUCH[®]

Dual mobility

Trapeziometacarpal prosthesis



Fully committed to **Hand & Wrist** only

www.kerimedical.com

TOUCH® Dual mobility trapeziometacarpal prosthesis


/ SPHERICAL CUPS

- 2 cup sizes (Ø 9 mm and 10 mm)
- 5 anti-rotation fins
- Spiked crown for stability
- Material / Stainless steel
- Dual coating / Porous titanium + HA



/ CONICAL CUPS

- 2 cup sizes (Ø 9 mm and 10 mm)
- Material / Stainless steel
- Dual coating / Porous titanium + HA



/ INSERT AND STRAIGHT NECKS

- Straight neck
- 3 heights / 6, 8 and 10 mm
- Head diameter 4 mm
- Insert material / Polyethylene
- Neck material / Stainless Steel



/ INSERT AND OFFSET NECKS

- 15° offset neck
- 3 heights / 6, 8 and 10 mm
- Head diameter 4 mm
- Insert material / Polyethylene
- Neck material / Stainless Steel



/ STEMS

- 5 stem sizes
- Anatomical trilobe section
- Material / Titanium
- Dual coating / Porous titanium + HA




THE + TOUCH®

- + Optimal Modularity

120 possible combinations to fit as close to the patient's anatomy

- + Make easier the revisions: each component can be revised independently

/ INSTRUMENT SET



KeriMedical Switzerland - HQ
 Route des Acacias, 45a
 1227 Geneva - Switzerland
 Tel: +41 (0) 58 255 01 30
 Fax: +41 (0) 58 255 01 40
 customerservice@kerimedical.com

KeriMedical France
 ActiPro - 219 rue Laszlo Biro
 74160 Archamps - France
 Tel: +33(0) 4 50 85 36 10
 Fax: +33 (0) 4 80 16 05 10
 adv@kerimedical.com

KeriMedical Germany GmbH
 Friedrich-Koenig-Str. 3
 55129 Mainz - Germany
 Tel: +49 (0) 6131 617 69 26
 Fax: +49 (0) 6131 617 69 28
 officegermany@kerimedical.com

