

TOUCH[®]

Dual mobility

Trapeziometacarpal prosthesis



TOUCH®

Dual mobility trapeziometacarpal prosthesis

/ SPHERICAL CUPS

- 2 cup sizes (Ø 9 mm and 10 mm)
- 5 anti-rotation fins
- Spiked crown for stability
- Material / Stainless steel
- Dual coating / Porous titanium + HA



CTO09

CTO10

/ CONICAL CUPS

- 2 cup sizes (Ø 9 mm and 10 mm)
- Material / Stainless steel
- Dual coating / Porous titanium + HA



CTO109

CTO110

/ INSERT AND STRAIGHT NECKS

- Straight neck
- 3 heights / 6, 8 and 10 mm
- Head diameter 4 mm
- Insert material / Polyethylene
- Neck material / Stainless Steel




NTO06

NTO08

NTO10

/ INSERT AND OFFSET NECKS

- 15° offset neck
- 3 heights / 6, 8 and 10 mm
- Head diameter 4 mm
- Insert material / Polyethylene
- Neck material / Stainless Steel



NTO156

NTO158

NTO1510

/ STEMS

- 6 stem sizes
- Anatomical trilobe section
- Material / Titanium
- Dual coating / Porous titanium + HA



STOX5

ST00

ST01

ST02

ST03

ST04

XS Stem

+ Optimal Modularity




144

possible combinations to fit as close to the patient's anatomy

- + **Make easier the revisions:** each component can be revised independently

THE + TOUCH®

/ INSTRUMENT SET



CE 1250

KeriMedical Switzerland - HQ

Route des Acacias, 45a
1227 Geneva - Switzerland
Tel: +41 (0) 58 255 01 30
Fax: +41 (0) 58 255 01 40
customerservice@kerimedical.com

KeriMedical France

ActiPro - 219 rue Laszlo Biro
74160 Archamps - France
Tel: +33 (0) 4 50 85 36 10
Fax: +33 (0) 4 80 16 05 10
adv@kerimedical.com

KeriMedical Germany GmbH

Friedrich-Koenig-Str. 3
55129 Mainz - Germany
Tel: +49 (0) 6131 617 69 26
Fax: +49 (0) 6131 617 69 28
officegermany@kerimedical.com

