

Nerve Regrowth and Protection tube made of CHITOSAN

REAXON[®]
DIRECT

PROTECTION TUBE

- Highly flexible
- Thin wall

- Made of polysaccharides
- Resorbable¹
- Transparent
- Easy to suture (hydrogel)

NERVE REGROWTH TUBE

- Flexible
- Does not collapse²
- Stable > 18 months³

REAXON[®]
NERVE GUIDE

Properties of CHITOSAN

- **Bioactive** - supports nerve regeneration equivalent to the autograft⁴
- **Antiadhesive** - inhibits scar tissue and neuroma formation⁵
- **Biocompatible** - does not induce any immune response⁶
- **Antibacterial** - limits infection risks⁷

REAXON® DIRECT - ULTRA FLEXIBLE PROTECTION TUBE

Reference	Inner diameter	Length
RD121	2.1 mm	14 mm
RD130	3.0 mm	14 mm
RD140	4.0 mm	14 mm
RD150	5.0 mm	14 mm
RD160	6.0 mm	14 mm

The diameter indicated is in a dry state. Once Reaxon® Direct has been soaked, its diameter increases up to 0.5 mm.

Reaxon® Direct is intended for use in repair of peripheral nerve discontinuities in patients with completely separated nerves that are treated by end-to-end coaptation.

REAXON® NERVE GUIDE - NERVE REGROWTH TUBE

Reference	Inner diameter	Length
RG321	2.1 mm	30 mm
RG330	3.0 mm	30 mm
RG340	4.0 mm	30 mm
RG350	5.0 mm	30 mm
RG360	6.0 mm	30 mm

The diameter indicated is in a dry state. Once Reaxon® Nerve Guide has been soaked, its diameter increases up to 0.5 mm.

Reaxon® Nerve Guide is intended for use in repair of peripheral nerve discontinuities in patients with completely separated nerves up to a nerve defect length of 26 mm.

¹ Internal testing at Medovent : Reports RDV-001 and RDV-002

² Internal testing at Medovent : Reports NG0064 and NG0090

³ Internal testing at Medovent : Report NG0086

⁴ Haastert-Talini et al. Chitosan tubes of varying degrees of acetylation for bridging peripheral nerve defects. Biomaterials 2013; 34: 9886-9904.

⁵ Marcol et al. Reduction of post-traumatic neuroma and epineural scar formation in rat sciatic nerve by application of micro-crystallic chitosan. Microsurgery 2011; 31: 642-649.

⁶ The biocompatibility of Reaxon® Nerve Guide has been proven in accordance with the ISO 10993 standards.

⁷ Rabea et al. Chitosan as antimicrobial agent: application and mode of action. Biomacromolecules 2003; 4: 1457-1465.

Note: Bench testing is not necessarily indicative of clinical performance.



Manufactured by Medovent, company of the KeriMedical Group

KeriMedical SA
Switzerland (HQ)
 Route des Acacias, 45a
 1227 Geneva - Switzerland
 Tél: +41 58 255 01 30
 Fax: +41 58 255 01 40
 customerservice@kerimedical.com

KeriMedical SAS
France
 ActiPro - 219 rue Laszlo Biro
 74160 Archamps - France
 Tél: +33(0)4 50 85 36 10
 Fax: +33 (0)4 80 16 05 10
 adv@kerimedical.com

KeriMedical GmbH
Germany/Austria
 Friedrich-Koenig-Str. 3
 55129 Mainz - Germany
 Tel: +49 (0) 6131 617 69 26
 Fax: +49 (0) 6131 617 69 28
 officegermany@kerimedical.com

KeriMedical
 www.kerimedical.com