



KeriFlex[®]

Silicone Finger Joint Implants

PIP - MCP




Fully committed to **Hand & Wrist** only

www.kerimedical.com


KeriFlex®

Silicone Finger Joint Implants PIP - MCP

Replaces joints affected by rheumatoid arthritis, osteoarthritis or traumatic arthritis



Test lab



Flexion / extension test of Implant

DESIGN & RESISTANCE VERIFIED

- Full flexion / extension (0° - 90°)
- Verified resistance up to 10 million cycles (test lab)
- Respect of joint angulation at rest : PIP 15° and MCP 30°

PIP 0

Size adapted to small morphologies

The new generation
of Finger Joint Implants

....

KeriFlex® implants provide
strength, flexibility
and durability

....

Medio-lateral stability
and joint mobility

....

High performance
medical silicone
for high mechanical
resistance
to physiological efforts

13 SIZES AVAILABLE

PIP - 6 sizes of implants

MCP - 7 sizes of implants


0°

■ Anatomic center of rotation is respected

■ Hinge design adapted to the joint, limiting conflicts with the extensor mechanism

■ Resistance to natural compression stress of the joint

90°



A = 9 mm
B = 6,5 mm

1 H = 6 mm

A = 10 mm
B = 7 mm



2 H = 6,9 mm

A = 10,8 mm
B = 7,8 mm

3 H = 7,6 mm

A = 12 mm
B = 8,5 mm

4 H = 8,2 mm



A = 10 mm
B = 9,4 mm

10 H = 7,5 mm

A = 10,8 mm
B = 10,2 mm

20 H = 8 mm

A = 12 mm
B = 11,4 mm


30 H = 9 mm

A = 13 mm
B = 12,4 mm

40 H = 10 mm

A = 13,8 mm
B = 13,5 mm

50 H = 10,5 mm



PIP & MCP INSTRUMENTATION

KeriFlex®


OPTIMISED & SIMPLIFIED
INSTRUMENTATION

Instruments


2 progressive rasps per indication




Set PIP



Set MCP



Set PIP + MCP



Validated sterilisation cycle 134°/18 min and 134°/3 min according to ANSI/AAMI ST79 and EN ISO 17665

ORDERING INFORMATION

PIP Implants

| Catalog | | |
|------------|----------------------------|----------|
| Reference | Designation | Pattern |
| 210-P20000 | PIP 0 Finger Joint Implant | ● Black |
| 210-P20001 | PIP 1 Finger Joint Implant | ● Orange |
| 210-P20002 | PIP 2 Finger Joint Implant | ● Green |
| 210-P20003 | PIP 3 Finger Joint Implant | ● Yellow |
| 210-P20004 | PIP 4 Finger Joint Implant | ● Red |
| 210-P20005 | PIP 5 Finger Joint Implant | ● Blue |

MCP Implants

| Catalog | | |
|------------|-----------------------------|----------|
| Reference | Designation | Pattern |
| 210-P10000 | MCP 00 Finger Joint Implant | ● Blue |
| 210-P10010 | MCP 10 Finger Joint Implant | ● Purple |
| 210-P10020 | MCP 20 Finger Joint Implant | ● Orange |
| 210-P10030 | MCP 30 Finger Joint Implant | ● Green |
| 210-P10040 | MCP 40 Finger Joint Implant | ● Yellow |
| 210-P10050 | MCP 50 Finger Joint Implant | ● Red |
| 210-P10060 | MCP 60 Finger Joint Implant | ● Black |

CE 1250

KERI MEDICAL SA

Route des Acacias, 45a
1227 Geneva - Switzerland
Tel: +41 58 255 01 30
Fax: +41 58 255 01 40
customerservice@kerimedical.com

KERI MEDICAL FRANCE

34 rue Antoine Redier
74160 Archamps - France
Tel: +33(0)4 50 85 36 10
Fax: +33 (0)4 80 16 05 10
adv@kerimedical.com

Distributed by:

KeriMedical